Summary Results of a Blinded Study on the Effectiveness and Efficiency of using SPERM HY-LITER™ to Screen Sexual Assault Evidence for Sperm.

Jennifer Old* Ph.D., Marisa Fahrner** MS, Jie Wu* Ph.D., Christian G. Westring*** Ph.D., P.W. Boonlayangoor* Ph.D. and Karl Reich* Ph.D. *Independent Forensics, **Commonwealth Biotechnologies, Inc., ***Department of Forensic Medicine, Section for Forensic Genetics

DNILM

Yes (cont)

Vaginal Introitus

Vaginal Fornix

No

No

No

Low Copy Samples (1-10 spermatozoa per prep)

Standard Samples, Dense Epithelia

No

Abstract

The identification of sperm is often a requirement of forensic DNA analysis for the processing and evaluation of sexual assault evidence. **SPERM HY-LITER™** is a new immunofluorescent detection kit specifically designed for forensic DNA laboratories. This kit uses a fluorescent tagged mouse monoclonal antibody with unique specificity for human sperm heads. The sensitivity of the test is such that a single sperm is easily and quickly identified on a microscope slide.

Forensic laboratories are naturally, and appropriately, skeptical of new and emerging technologies. At the request of the Department of Forensic Genetics of the Danish National Institute for Legal Medicine (the centralized forensic DNA laboratory of Denmark), Independent Forensics agreed to perform a blinded study on the specificity, sensitivity and work flow efficiency of **SPERM HY-LITER™** for the microscopic screening of sperm from sexual assault evidence.

The Danish National Institute for Legal Medicine (DNILM) prepared a duplicate series of swabs and stains with forty-four (44) samples prepared from a variety of sources including duplicate cuttings from sexual assault evidence kits retained at the Danish National Institute for Legal Medicine. All samples were blinded and number coded to both laboratories. An identical protocol for the extraction, processing and staining of all samples was chosen prior to testing – the provided SPERM **HY-LITER™** protocol was followed with the addition of a short, waterbath sonication extraction step and Spin-Eze™ aided recovery of the resultant cell pellet. The Department of Forensic Genetics used both SPERM HY-LITER™ staining and their standard histological staining (i.e., H&E) for the evaluation of all samples.

Upon completion of the testing and the exchange of data, sample codes were revealed. Sample types tested and evaluated included negative controls, low sperm count samples, dense samples, epithelial, vaginal, washed semen stains, mixtures of blood, menstrual blood, soil, saliva, and bacterial with seminal fluid all formed part of the tester set.

Data comparison of the three sets of samples (SPERM HY-LITER™ staining at IFI, **SPERM HY-LITER™** staining at DNILM, and H&E staining at DNILM) revealed that **SPERM HY-LITER™** was clearly the most sensitive of the staining methods.

The most dramatic result of the study is the time savings seen with **SPERM HY-LITER™:** the two weeks required by DNILM to stain and scan the forty-four samples using their standard H&E histological method was reduced to *three (3) days* with **SPERM HY-LITER™.** This represents an 80% time saving with increased sensitivity as compared to standard forensic technique.

Methods

Extraction

All samples were extracted using an identical protocol with standardized volumes, times and conditions. Briefly, swabs and cuttings were soaked in phosphate buffered saline (PBS) for 1-2 hrs at room temperature, and incubated in a water bath sonicator for twenty (20) minutes and the extract recovered, and cells pelleted using Spineze and centrifugation. Cell pellets were resuspended in 50-200 µL of PBS.

All samples were stained using **SPERM** standard DNILM method (H&E).

Documentation

Slides were scored visually and communicated between IFI and DNILM before coded information was revealed.

Total time –

Negative Controls

Buccal Swabs

Vaginal Fornix

Rectum

Epithelial

Vaginal Introitus

Vaginal with blood

Buccal Swabs

Vaginal Fornix

Rectum

Epithelial

Buccal Swabs

Vaginal Introitus

Vaginal with Blood

Vaginal Fornix

Rectum

Epithelial

Vaginal Introitus

Vaginal with blood

Staining


HY-LITER™ as described by the manufacturer (IFI) and using the

Time for Analysis

IFI: three (3) days DNILM: two (2) weeks

Time per sample – IFI: 32 min DNILM: 109 min

Results



Conclusions

Negative Controls

IFI - All negative controls correct DNILM – one false positive

Low Copy Samples

IFI - 1 of 6 **DNILM** – 1 of 6

Standard Samples

IFI - 5 of 6 **DNILM** – 3 of 6

Degraded Samples

IFI- 4 or 6 **DNILM** – 4 of 6

Experimental Samples

IFI-14 of 20 **DNILM 12 of 20**

Sensitivity:

SPERM HY-LITER™: 30 of 44 **DNILM: 25 of 44**

Specificity:

SPERM HY-LITER™: 0 false positives DNILM: 1 false positive

Efficiency:

SPERM HY-LITER™: 80% time saving Three (3) times faster screening of

